



The Medical Power of Light

ClearSteps™ Laser Onychomycosis Treatment

- Complete fungal elimination
- Non-invasive
- No chemicals or oral medication
- Stimulates the natural growth and immune processes of the body
- Extremely fast, effective and easy to perform
- Safe and patient-friendly treatment

Fotona
choose perfection

www.fotona.com


89805 CE ENG/2

1




Step 1:
File down the over-grown nail in order to fully penetrate into the treatment site.

2



Step 2:
Pass the light beam over the nail to eliminate the fungi.

3



Step 3:
Treat the skin tissue surrounding the nail to prevent recurrence.

What is ClearSteps™ Onychomycosis Treatment?

ClearSteps is a revolutionary method for treating onychomycosis in a patient-friendly, quick and effective way with Nd:YAG laser light. Laser light heats evenly throughout the depth of effected nail and skin tissue, effectively weakening and killing parasitic fungi which have infected the patient's nail.

How does ClearSteps™ work?

Deep, pulse-profile heating of the nail bed stimulates the killing of parasitic fungus. The natural growth and immune processes of the body are then able to restore the nail to its pristine state. Compared to traditional methods the effects of laser light spread evenly throughout the tissue and are not subject to either the limits of chemical diffusion, as are topicals, or the induction of hepatotoxic side effects, as is the case with oral medications.


The procedure can be performed in **three simple steps**:

Step 1: File down the Nail. During onychomycosis the nail will become deformed and overgrown, it may also thicken. Filing down the nail allows laser light to fully penetrate into the treatment site.

Step 2: Administer Nd:YAG laser pulses in a circular pattern covering the entire nail. Slightly overlap the beam so that the nail is fully covered. The patient should feel a slight heating sensation throughout the nail.

Step 3: Target Nd:YAG laser pulses around the perimeter of the nail bed and around the finger or toe. Often parasitic fungus has worked its way beyond the nail bed and needs to be eliminated in all the tissue surrounding the nail to prevent recurrence.

The treatment is usually performed four times at one week intervals. All of the nails should be treated, during every treatment, just in case. Laser light can also effectively treat fungus which has infected skin tissue.



The Nd:YAG laser penetrates the nail efficiently, safely and evenly heating the nail bed to an appropriate temperature which results in the elimination of the invading fungi and the restoration of the nail.



In just 12 months this patient was completely cleared from severe onychomycosis with just the use of Fotona's laser.

Why is the Dynamis Line perfect for ClearSteps™?

The pulse profile, speed and power of lasers in the Dynamis Line allow a patient's nails to be treated quickly and effectively. Unparalleled speed makes the procedure convenient and

comfortable for the patient. In 15 minutes your patients can be on the way to enjoying healthy, aesthetically pleasing nails without the dangerous side-effects associated with anti-fungal drugs.

Exclusive advantages of ClearSteps™ for you and your patients

ClearSteps Onychomycosis Treatment is the most effective and least invasive treatment of fungal nail infection. The treatment is extremely fast and easy to perform and has an unsurpassed treatment success rate compared to traditional methods.

An important advantage of the ClearSteps treatment method is that no oral medications or chemicals are used. Therefore, there are no dangerous side effects as there are with conventional treatments. ClearSteps is the most patient-friendly and non-invasive treatment of onychomycosis.

Getting Started with ClearSteps™

It's easy to begin offering ClearSteps Onychomycosis Treatment to your patients. Training is provided through Fotona's partnership with the Laser and Health Academy (www.laserandhealth.com). It is available stand-alone or as part of one of the many comprehensive training programs provided by LAHA. After training you will be able to confidently treat fungal infections the natural way.

To learn more about **ClearSteps™** and what the **Dynamis Line** can do for your practice contact **Fotona** at info@fotona.com today.